

# 200 SERIES GC & AUTOSAMPLER

Mult-Level Temperature Programming  
Electronic Carrier Gas Control  
Split/Splitless Injector  
Choice of FID, TCD and ECD Detectors  
Capillary Column Designed Instrument  
Packed Column Adaptor available  
Small foot Print 41cm x 34cm  
Lightweight Only 7.5Kg

22 Position Autosampler  
2 Wash Bottles  
Method Settable by dip-switches  
Self Aligning carousel



**CSI**

**Advanced  
Chromatography  
Systems**